

# Exploration of Patent Chemistry by Fuzzy MCS-led Fragment Decomposition

*Richard Sherhod*

[r.sherhod@vernalis.com](mailto:r.sherhod@vernalis.com)



## • Expertise


- Fragments and structure-based drug discovery  
(Protein Science, Structural Biology, Chemistry)

## • Therapeutic areas

- Oncology, CNS, infectious diseases

## • Location

- Based in Granta Park, outside Cambridge, UK



- **Our aims**
  - Our solution
- **Patent data-mining**
  - SAR table extraction
  - Data curation
- **Fuzzy MCS analysis**
  - Maximum common substructure mining
  - Fuzzy MCS
  - Recursive MCS mining
  - Scaffold-tree formation
- **Fragment decomposition**
  - MCS-led decomposition
  - Implementation
- **Conclusion**

## *Our Aims*



## The situation:

- New project with existing target-related patents
- ...or existing project with newly-published target-related patents
- We need to understand the chemical space covered by the patents

## Our aims:

- To summarise the coverage of exemplified structures
- To extract and present relevant SAR data


- KNIME based application for exploring patent chemistry
  - Accessed via the KNIME Web Portal
- KNIME workflows for preparing data and presenting results
  - Patent processing workflows (admin):
    - Structures and data extraction and curation
    - Recursive fuzzy MCS mining
  - Interactive results workflow
    - MCS tree formation for visualisation
    - Fragment decomposition for chosen MCS



## *Patent Data Mining*



- Patent data mined with NextMove Software's LeadMine
  - Processes HTML, XML, raw text etc.
  - Automatic structure extraction
    - Text-to-structure – IUPAC, generic names, abbreviations
    - CDX-to-structure – Ambiguities, drawing errors e.g. floating alkanes




May, J., Lowe, D. & Sayle, R., 2016. Sketchy Sketches: Hiding Chemistry in Plain Sight. Seventh Joint Sheffield Conference on Chemoinformatics. Available at: <http://cisrg.shef.ac.uk/shef2016/talks/poster21.pdf> [Accessed September 14, 2017].

- Patents accessed and processed by PatFetch web service
  - Patent archives stored locally



- Structures extracted from:
  - Names
  - CDX images
  - R-group tables
- Table data – activity, properties etc.
  - New feature of LeadMine
- IDs and data associated with structures



**3-methyl-N-(4'-((1-methylpiperidin-4-yl)oxy)-[1,1'-biphenyl]-3-yl)benzamide (6d)**



**TABLE 2**  
Activity of Para-Meta Biphenyl Core with Substituted Benzylamide

Entry	R	SKBr3 (IC <sub>50</sub> , μM)	MCF-7 (IC <sub>50</sub> , μM)
6b	H	18.86 ± 0.95	12.02 ± 0.57
6c	p-CH <sub>3</sub>	5.27 ± 0.29 <sup>a</sup>	3.92 ± 0.13
<b>6d</b>	m-CH <sub>3</sub>	11.38 ± 1.37	7.73 ± 1.90
6e	p-t-butyl	1.51 ± 0.31	3.45 ± 0.02
6f	p-methoxy	10.1 ± 0.93	5.52 ± 0.01
6g	m-methoxy	8.36 ± 1.35	4.50 ± 0.46
6h	p-Cl	3.63 ± 1.03	2.23 ± 0.05
6i	m-Cl	4.29 ± 0.43	2.11 ± 0.42
6k	o-Cl	7.87 ± 0.48	5.17 ± 0.49
6l	p-Br	1.94 ± 0.11	0.88 ± 0.07
6m	3,4-dichloro	2.24 ± 0.11	2.17 ± 0.37
6n	2,4-dichloro	5.91 ± 0.15	3.93 ± 0.47
6o	3,5-dichloro	4.23 ± 0.09	3.72 ± 0.15
6q	-(2-naphthoyl)	2.09 ± 0.34	1.66 ± 0.27
6p	-(1-naphthoyl)	1.64 ± 0.13	1.10 ± 0.17

US-20160272584-A1

- Extracting clean patent data is difficult
  - Patent tables can contain errors:
    - Typos
    - OCR errors
    - Missing values
    - Inconsistent labelling schemes
  - Structure names and CDX images can disagree

4-{1-[3-(4-Bromo-benzoyl-amino)-phenyl]-  
3-pyrrolidin-1-yl-propyl-amino}-  
quinazoline-8-carboxylic acid amide



US 8637532 C00342

- No two patents are alike


Lowe, D., Senger, S. & Sayle, Ro., 2017. Automatic extraction of bioactivity data from patents. *253rd ACS National Meeting, San Francisco, CA, USA*. Available at: <https://www.slideshare.net/NextMoveSoftware/automatic-extraction-of-bioactivity-data-from-patents-74402139> [Accessed September 14, 2017].

- Fully automated extraction is unattainable (for now)
  - Structures and data are interactively assessed
  - Data may be exported for manual editing
  - Modified data is semi-automatically validated


## *Fuzzy MCS Analysis*



- The largest substructure common to a set of structures




- ...or disconnected substructures



- Many implementations and uses in cheminformatics
- Traditionally only exact atom and bond matches allowed


- MCS with variation in atoms and/or bonds




- Algorithm available in RDKit
  - Contributed by Andrew Dalke in 2012
  - RDKit MCS KNIME node



- Set of fuzzy MCS generated from patent structures
  - MCS generated for a range of coverage thresholds
    - i.e. MCS that represent 10%, 20%, 30%... 100% of input structures
  - Structures not covered by MCS are used to generate new MCS




- Resulting collection of MCS arranged into trees
- Hierarchical relationship between MCS coverage
- Fuzzy (overlapping) hierarchical clustering






- Custom code used to generate and visualise tree structure
  - Tree represented as JSON object
- Tree presented as interactive view in KNIME Web Portal
  - Tree visualised in D3.js
  - Crude POC




# *Fragment Decomposition*




- Break structures down into categorised fragments
  - R-groups
  - Chains – terminal or linkers
  - Rings – fused systems
- R-group decomposition
  - No scaffold
  - Suited to FBDD




- Fuzzy MCS used to define fragment framework
- Fuzzy MCS decomposed into:
  - R-groups, Rings, Chains
  - Rings with variable features
  - Chains with variable features
- Fuzzy MCS fragments used to decompose structures



- In-house algorithm developed with RDKit Java API
- Fragment Decomposition KNIME node released internally
  - Thanks Steve Roughley!





- Fuzzy MCS is fragmented
  - Break all non-aromatic bonds between ring and non-ring atoms
    - SMARTS: [!R0]!@&!:




- Fuzzy MCS is fragmented
  - Break all non-aromatic bonds between ring and non-ring atoms
    - SMARTS: [!R0]!@&!:
- Fragments are categorised
  - R-group, Ring, Chain, qRing, qChain
- Fragments given canonical identifiers







- Fragments are mapped to input structures




- Fragments are mapped to input structures




- Fused aromatic ring systems are expanded
  - Aromatic ring bonds are not broken
- Rings and chains are fragmented further
  - IDs of additional R-groups are canonicalised in a later process




- qRings and qChains are either:




- qRings and qChains are either:
  - Expanded to incorporate additional features



- qRings and qChains are either:
  - Expanded to incorporate additional features




- qRings and qChains are either:
  - Expanded to incorporate additional features
  - Fragmented further into new R-groups





- qRings and qChains are either:
  - Expanded to incorporate additional features
  - Fragmented further into new R-groups



# Fragment Decomposition






















## Implementation – Visualisation

### Results

Please select results for export to JChem for Excel


Show  entries

Search:

<input type="checkbox"/>	ID	Structure	MCS ID	MCS	qRing1	qRing2	qRing3	Ring1	qChain1
<input type="checkbox"/>	19a		MCS-3						
<input type="checkbox"/>	19b		MCS-3						
<input type="checkbox"/>	19c		MCS-3						

# *Demonstration*



 Open for Innovation **KNIME** gadolinium

AdministrationSettingsLogout

Explore patent chemistry

☐ Mail notification on completion

Start ▶

©2017 KNIME.com AG - Version 4.5.0

## *Conclusion*



## Developed tools/algorithms for:

- Extracting structures and data from the patents
- Summarising exemplified structures as fuzzy MCS
- Showing hierarchical relationship between fuzzy MCS
- Performing fragment decomposition driven by fuzzy MCS
- Presenting results to users

## Future work:

- More automation of patent document processing
  - Learn lessons from processing more patents
- Multi-parametric SAR/SPR analysis
  - Process and compare multiple tables, e.g. binding and stability data tables
- Integration with other services/tools

*Thank you!*



THE QUEEN'S AWARDS  
FOR ENTERPRISE:  
INTERNATIONAL TRADE  
2014